SIGNATURE SERIES



Regal PR751TY

The Regal is a feature-rich lift/recliner designed for people who like to keep things within hands' reach. Its dual flip-open arms conceal extra storage space, a cup holder, and hideaway tray. A luxurious three-pillow waterfall back can be adjusted, providing more comfort options. Stylish and superbly constructed, the Regal will allow you to kick back and relax in first-class comfort.





Regal



Pogal

Madel Number	PR-751TY
Model Number	, 5
Size	Medium
Weight Capacity	375 Lb.
Recline Positions	3
Standard Fabrics	5
Fabric Options	ALL
Chaise Seat	YES
Coil Springs	YES
Back Type	3 Pillow
Head Rest Cover	YES
Arm Covers	NO
Battery Back Up	YES
Overall Width	36"
Overall Height	44"
Width Between Arms	22"
Floor to Top of Seat	20"
Seat Depth	21.5"
Seat to Top of Back	28"
Seat to Top of Arm	8"
Seat Angle Raised	18°
Floor to Top of Foot Rest Reclined	24"
Floor to Top of Back Reclined	33"
Back Angle Closed	110°
Back Angle Reclined	140°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	180°
Distance Required from Wall Reclined	18"
Extension Reclined from Back of Seat to End of Foot Rest	42"
Life Time Warranty, Chair Frame,	YES
Lift Frame, Recline Mechanism	
Warranty on Electrical Parts	3 Yrs.
Pro Rated Years 4 Thru 7	YES
Warranty on Mechanical Labor	3 Yrs.
Warranty on Electrical Labor	1 Yr.
Removable Back Standard	YES
Smart Tek System Standard	YES
Customizable	NO
Left Hand Control Option Available	YES
Extra Pocket Option Available	YES
Heat and Massage Option Available	YES
Foot Rest Extension Option Available	NO
Straight Lift Option Available	YES
User Height Suggested	5'4"- 5'10"

Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



Battery Backup System

Battery Backup System

Batteries not included.
Battery Backup System is
designed to bring the chair
from the reclined to the
seated position in the event
of a power outage.



Smart Tek

In the event of a malfunction, our patented SmartTek™ diagnostic system enables you to get your lift chair up and running quickly.



Golden Craftsmanship

Each luxury lift & recline chair is hand crafted in our facility using decades of motion furniture experience.



Added Comfort and Safety

Designed to perform year after year using advanced technology backed by the best warranty.



Smart Engineering

Smooth, quiet lifting system offers a secure lift for the best performance.



Literature is current at time of printing Golden Technologies reserves the right to make changes to the product or literature at any time.





