

New and Improved Controller!



Seat Colors  
Black Vinyl

Body Colors  
Red  
Blue

# golden COMPASS<sup>TM</sup> SPORT

GP605

The Golden Compass Sport features the new Dynamic Linx controller! The new Linx offers a better ride experience using patented Dynamic Load Compensation technology. This new technology improves the terrain-climbing ability of the Compass Sport's patented springless articulating chassis. Learning to drive the Compass Sport is easier and a more natural experience than before!



## Golden Compass™ Sport Specifications

<b>Model</b>	<b>GP605</b>
Model Name	Compass Sport
Medicare Code	K0823
Maximum Incline	6°
Weight Capacity	300 lbs
Drive Wheels	Center
Maximum Speed	4 mph
Operating Range <sup>1</sup>	24 miles
Ground Clearance	3.25"
Turning Radius	19.5"
Type Batteries	2-U1
Colors	Red, Blue
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	40 Amp
<b>Dimensions:</b>	
Length	39"
Width (outside of tires)	24"
Height (ground to top of back of seat)	38" - 40"
Height (ground to top of back of headrest)	45" - 50.5"
Ground to Top of Seat	21" - 23"
Number of Seat Height Positions	3
Size of Increments	1"
<b>Weight of Unit: (assembled)</b>	
Base	92 lbs
Seat with Arms	46 lbs
Batteries	22 lbs
<b>Tires:</b>	
Center	10"
Casters	6"
<b>Standard Seat:</b>	
Back Height - no headrest	19"
Back Height - with headrest	26.5" - 30"
Width x Depth	18" x 18" / 20"
Color	Black Vinyl
<b>Optional Seats:</b>	
Back Height - no headrest	19"
Back Height - with headrest	26.5" - 30"
Width x Depth	16", 18", 20", or 22" x 16", 18" or 20"
Color	Black Vinyl
<b>Footrest:</b>	
Height adjustable	3" - 7"
Angle adjustable	0 - 35 degrees
<b>Warranty:</b>	
Frame	Lifetime
Drive Train	13 Months
Electronics <sup>2</sup>	13 Months
Batteries	12 Months

<sup>1</sup> Battery range at 200 lbs. and will vary due to rider weight, drive surface, terrain & battery type.

<sup>2</sup> Electronics warranty excludes batteries.

Armrests are height, width and angle adjustable



40 amp; fully programmable controller

Padded flip-up arm support

Exclusive fold-flat seat hinge

Center wheel drive

Foot rest adjustable:  
-height  
-length  
-angle

## 2 gorgeous colors



Red



Blue

## Available Accessories:

- Crutch Holder Tube
- Cane Holder Tube
- Oxygen Tank Holder
- Quad Cane Holder
- Walker Holder
- Swing-away Desk Arm
- Swing-away Footrests
- Stump Support
- Elevating Leg Rests
- Power Chair Cover
- Rear Basket
- Pack N' Go
- Van/pan Seat
- EZ Lock System
- Universal Cup Holder
- Armrest Bag